



VitaTaq® DNA Polymerase is a robust and reliable Taq DNA polymerase suitable for all common PCR applications like colony PCR, cloning applications, high-throughput PCR and routine PCR.

PRODUCT	SIZE	SKU
VitaTaq®DNA Polymerase + 10X Buffer	500 units	PCCSKU1001
	1000 units	PCCSKU1002
	5000 units	PCCSKU1003
2X HS Mastermix Clear	2 ml	PCCSKU1004
	5 ml	PCCSKU1005
	10 ml	PCCSKU1006
2X Mastermix Gold	2 ml	PCCSKU1007
	5 ml	PCCSKU1008
	10 ml	PCCSKU1009

#### STORAGE CONDITIONS

Store all components at -20°C and avoid repeated freeze and thaw cycles.

#### ADDITIONAL MATERIALS REQUIRED

Nuclease free dH<sub>2</sub>O

Nuclease free PCR tubes / plates & sealing options

Thermocycler

PCR Primer (10 µM each)

dNTP Mix (10 mM each)<sup>1</sup>

template DNA

<sup>1</sup> (not required for 2X Mastermix Clear/Gold)

#### REACTION SETUP

1) Thaw all reaction components completely and mix gently to ensure even distribution of all components. Prepare the reaction on ice in a sterile, nuclease free tube and mix gently after addition of the polymerase. Collect all liquid at the bottom of the tube by a quick spin.

COMPONENT	VOLUME	FINAL CONCENTRATION
10X PCR Buffer	5 µl	1X
dNTP Mix (10 mM each)	0.5 µl	0.1 mM each
Primer 1 (10 µM)	1 µl	0.1 µM – 0.5 µM
Primer 2 (10 µM)	1 µl	0.1 µM – 0.5 µM
VitaTaq® DNA Polymerase (5 U/µL)	0.2 µl	1 U
template DNA	1 µl	<1 µg
dH <sub>2</sub> O	to 50 µl	

2) Keep the reactions on ice until transfer to the thermocycler, then cycle according to these guidelines:

STEP	CYCLES	TEMPERATURE	DURATION
Initial Denaturation	1	94°C	5 minutes
Amplification	30-35	94°C	30 seconds
		$T_m - 5^\circ\text{C}$	30 seconds
		72°C	1 minute / kb
Final Extension	1	72°C	5 minutes
Hold	1	4°C	

3) Analyse the amplification reaction by gel electrophoresis using an acrylamide or agarose gel of appropriate percentage.

## 2X HS Mastermix Clear / Gold

For routine applications, VitaTaq® DNA Polymerase is available as a 2x mastermix formulation; all there is to add is your primers and template DNA.

2X HS Mastermix Gold furthermore contains a tracking dye that enables direct gel loading without addition of a separate loading dye. For increased fidelity, 2X Mastermix Clear/Gold comes with a HotStart feature that allows reaction setup at room temperature.

### REACTION SETUP

1) Thaw all reaction components completely and mix gently to ensure even distribution of all components. Prepare the reaction in a sterile, nuclease free tube and mix gently. Collect all liquid at the bottom of the tube by a quick spin.

COMPONENT	VOLUME	FINAL CONCENTRATION
2X HS Mastermix Clear/Gold	25 µl	1X
Primer 1 (10 µM)	1 µl	0.1 µM – 0.5 µM
Primer 2 (10 µM)	1 µl	0.1 µM – 0.5 µM
template DNA	1 µl	<1 µg
dH <sub>2</sub> O	to 50 µl	

2) Transfer the reactions to the thermocycler, then cycle according to these guidelines:

STEP	CYCLES	TEMPERATURE	DURATION
Initial Denaturation	1	94°C	5 minutes
Amplification	30-35	94°C	30 seconds
		$T_m - 5^\circ\text{C}$	30 seconds
		72°C	1 minute / kb
Final Extension	1	72°C	5 minutes
Hold	1	4°C	

3) Analyse the amplification reaction by gel electrophoresis using an acrylamide or agarose gel of appropriate percentage.